



## Tellus 3+

### Tellus 3+

The Tellus 3+ is a Fujitsu tablet computer that has been adapted to create a very robust, wheelchair mountable dynamic screen communication aid.

#### What does it include?

- 512MB SDRAM (upgradeable to 4GB), Intel® Core™ Duo Processor Ultra Low Voltage U2500 (1.20 GHz, 2MB L2 cache, 533 MHz FSB);
- 60 GB Hard Drive
- Windows XP™ Tablet Edition
- Integrated amplified speakers
- Integrated wireless IEEE802.11a/b/g, 2 x USB 2.0, Type II PC Card, IEEE 1384 Firewire, Bluetooth v2.0
- 5 Programmable buttons on the front of the Tellus 3+ - can be used for volume control etc
- DaeSSy Compatible Mounting Plate
- MindExpress & Eurovocs suite with extra desktop licence
- USB DVD/CDRW Drive, Infra-red keyboard, USB mouse, USB Memory Stick and Soft Carry Case
- 3 Years Warranty and 3 Years Screen Damage cover
- 1 Day of Familiarisation Training



Colour: Blue

Tellus 3+ - Colour options below



#### Options Include...

Finger Touch Screen (resistive)...True finger touch access  
GEWA Module for environmental control  
External battery or wheelchair power cable  
Switch or Joystick interface

#### Software...

Tellus 3+ is the perfect platform to run  
\* [Mind Express](#) symbol communication software  
\* [Eurovocs Suite](#) for text-based communication  
\* a range of Windows based software programs

#### PC based...

Allows videos, photos, drawings, music and sounds to be imported and used within the system.

#### Built in amplified speakers...

Good sound quality and audibility. High quality voice of RealSpeak™ also available for use with MindExpress™ and Eurovocs Suite™.



Tellus 3+ on DaeSSy Deskmount

#### Bluetooth...

Using the integrated Bluetooth module, a suitable mobile phone and the Mind Express™ SMS module, the Tellus 3+ can be used to send and receive text messages.

Bluetooth can also be used to enable wireless printing from the Tellus 3+.

#### Battery Life...

The internal battery lasts for 3½ hours. This can be extended by the use of an External battery. This will give 8 hours continuous use including the 3½ hours of the Tellus 3+ battery. The battery can be slotted into a special mounting plate on the back of the Tellus 3+. Instead of the external battery an optional wheelchair power cable can be used to run the Tellus 3+ from the batteries in an electric wheelchair.